



Revolutionizing Multiplex PCR: A Proof-of-Concept Study Demonstrating Four-Target Detection (Influenza A and B, RSV, and IC) with Two Fluorescence Channels in One Tube

JungKwon Yoon¹, Jaehyeon Lee²

¹ Amplidea, Jeonju, Republic of Korea

² Department of Laboratory Medicine, Jeonbuk National University Medical School and Hospital, Jeonju, Republic of Korea

Backgrounds

- ❖ Multiplex real-time PCR is a cornerstone of respiratory-virus diagnostics.
- ❖ Limited number of fluorescence channels constrains simultaneous target detection.
- ❖ **FluoroSplit** technology: signal-splitting platform that resolves two targets within a single channel (figure 1)
 - On a five-channel real-time PCR platform, up to 10 distinct targets can be simultaneously detected using the FluoroSplit
 - Expands multiplex capacity without additional hardware.
- ❖ This is a Proof-of-concept study
 - four-target assay detecting Influenza A (FluA), Influenza B (FluB), RSV, and an internal control (IC) in one tube.
 - Presenting proof-of-concept data and a preliminary clinical evaluation of prototype reagent.

Methods

- ❖ A four-target multiplex assay using FluoroSplit was developed.
 - ❖ Platform: CFX96 Touch instrument (Bio-Rad, USA).
 - ❖ Channels:
 - ❖ FAM → FluA and FluB
 - ❖ Cy5 → RSV and IC
 - ❖ Acquisition: dual fluorescence acquisition at 62 °C and 84 °C.
- ❖ Specimens: 34 RNA extracts from clinical nasopharyngeal swabs (previously tested by Allplex Respiratory Panel, Seegene, Republic of Korea).
- ❖ Contrived mixtures (n=10): to confirm clear target separation within a single channel
 - ❖ 2 mixtures with FluA + FluB
 - ❖ 4 mixtures with FluB + RSV
 - ❖ 4 mixtures with FluA + RSV
- ❖ Statistical Analysis: Agreement statistics using linear-weighted kappa calculated with MedCalc 23.3.2.




Figure 1. This technology separates fluorescence signals by temperature, enabling each target to be detected independently in real-time PCR and producing separate Ct values. Unlike melting curve analysis, it requires no additional post-run analysis, allowing amplification of each target to be monitored in real time.

Results

- ❖ Concordance: Weighted kappa = 0.81 (95 % CI, 0.64-0.95).
 - Missed detections: 5 high-Ct specimens (3 FluA, 2 RSV).
- ❖ Contrived samples: 90 % differentiation achieved.
 - FluA and FluB separated in FAM.
 - RSV and IC separated in Cy5.
 - Internal control: Recovered in every reaction, confirming assay integrity.




Figure 3. Contrived mixture: FluA + RSV. The internal control amplifies in Quasar 670 at 62 °C (left). At 84 °C, distinct curves for FluA (FAM) and RSV (Quasar 670) are visible in the right.




Figure 2. Contrived mixture: FluB + FluA. In the left, FluB is observed in FAM at 62 °C and the IC in Quasar 670. The right figure shows the Influenza A (FluA) curve in FAM at 84 °C. Note how the FluB signal at 62 °C and the FluA signal at 84 °C remain cleanly separated without interference.

Conclusion

- ❖ FluoroSplit assay enables accurate detection of four respiratory targets using only two fluorescence channels.
- ❖ **Proof-of-concept: two targets can be distinctly resolved within a single channel.** On five-channel platforms (e.g., CFX96), approach could support decaplex assays.
- ❖ Operates via temperature-dependent signal acquisition only:
 - No specialized hardware required.
 - No post-run mathematical deconvolution required.
 - Can be implemented immediately on standard platforms (e.g., CFX96).
- ❖ Particularly advantageous for: 1) Resource-limited settings, 2) Syndromic panels constrained by instrument capability.
- ❖ Findings underscore potential of FluoroSplit to deliver high-efficiency molecular diagnostics with minimal fluorescence bandwidth requirements for clinical laboratory practice.